




# EntericBio



## REVOLUTIONISING MOLECULAR ENTERIC TESTING

CE IVD

No DNA/RNA  
extraction

No manual  
pipetting

Rapid results  
within 3 hours


High  
throughput



## SEROSEP




# WORKFLOW



\* Workflow for the EntericBio Dx Panel processes 30 patient samples and a +/- control

## 4 STEP PROCESS



- 1-46 patient samples processed per run (including +/- controls)
- Lyophilised PCR reagents
- Less than 30 mins hands on time
- Direct from stool sample
- Time to result: under 3 hours
- High throughput up to 240 samples in an 8 hour working day

# PANELS

### EntericBio® Dx

- Bacteria**
  - Salmonella
  - STEC
  - Shigella
  - Campylobacter
  - Yersinia
  - Vibrio
- Parasite**
  - Cryptosporidium
  - Giardia
  - Entamoeba histolytica

### Viral

- Viral**
  - Adenovirus
  - Astrovirus
  - Norovirus GI
  - Norovirus GII
  - Rotavirus
  - Sapovirus

### Clostridioides difficile

- Bacteria**
  - Toxin B (tcdB)

### CPE\*

- CPE**
  - OXA - 48
  - KPC
  - NDM
  - NIM
  - IMP
  - GES

\* Specimen type: Stool samples or rectal swabs (in liquid medium) and bacterial isolates

Annacotty Business Park, Annacotty, Co. Limerick, V94 FF83, Ireland  
+353 61 358 190 | info@serosep.com | www.serosep.com



Scan for further information

Distributed by: **EUROIMMUN**

EUROIMMUN Medizinische Labordiagnostika AG  
Seekamp 31 · 23560 Lübeck (Germany)  
Phone +49 451 2032-0 | Fax +49 451 2032-100  
www.euroimmun.com