


EURORealTime HSV-1/2

Quantitative direct detection of herpes simplex virus types 1 and 2 by real-time PCR

- For the detection of viral DNA in smears, CSF and EDTA plasma
- Simple quantification of the viral load in sample material based on ready-for-use standards
- Automated evaluation and documentation using the EURORealTime Analysis software

EURORealTime HSV-1/2


(Order number: MP 2530-####)




^{*} Test is not available in the USA

EURORealTime Analysis software

- Direct quantification of the viral load
- Reusable standard curves
- Simplified workflow through software support with automated generation of pipetting schemes
- Efficient process from standardised raw data analysis to documentation of results
- Compatible with different real-time PCR cyclers and LIS



(Order number: YG 0661-0101)



Analytical performance

- Lower detection limit (LoD): 3cp/µl (DNA copies per µl eluate) for both HSV subtypes
- No cross-reactions detected with DNA from various microorganisms commonly found in anogenital areas or closely related to HSV-1 or HSV-2